FT3

Free Triiodothyronine (CLIA)

Order Information

Catalog No.	Package Size
FT3111	2×50 tests
FT3112	2×100 tests

Intended Use

The CL-series FT3 assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of free triiodothyronine (FT3) in human serum.

Summary

Thyroxine (T4) and Triiodothyronine (T3) are thyroid hormones, which regulate normal growth and development by maintaining body temperature and stimulating calorigenesis.¹ The hypothalamic-pituitary-thyroid axis controls thyroid hormone synthesis, release, and action. Thyrotropin-releasing hormone (TRH) secreted from the hypothalamus stimulates the synthesis and release of thyrotropin or thyroid-stimulating hormone (hTSH). hTSH, in turn, stimulates the synthesis, storage, secretion, and metabolism of T4 and T3.

T3 is the major biologically active thyroid hormone, with a molecular weight of 651 daltons. Of the circulating T3, about 80% is formed from peripheral deiodination of thyroxine and 20% is secreted directly from the thyroid gland.²

In circulation, only 0.2%-0.4% of the circulatory total T3 is in equilibria as free, and others are bounded to thyroxine binding globulin (TBG), thyroxine binding pre-Albumin (TPBA) and albumin.^{3,4} Only free T3 can function.^{5,6}

Free T3 level is obviously elevated in hyperthyroidism, and mostly with the elevation of free T4.^{1,7,8} However, free T3 is elevated alone (T3 toxicosis) in about 5% of hyperthyroids.⁹ The determination of free T3 has the advantage of being independent of changes in the concentration and binding properties of the binding proteins. Thus, free T3 assay is always to be used as an aid in the assessment of thyroid status.

Assay Principle

The CL-series FT3 assay is a competitive binding immunoassay to determine the level of free T3.

In the first step, sample, monoclonal anti-T3 antibody -alkaline phosphatase conjugate are added into a reaction vessel.


After incubation, free T3 present in the sample binds to anti-T3 antibody. Then biotinylated T3 and microparticles coated with streptavidin are added to the mixture. Biotinylated T3 competes with free T3 in sample for binding site of anti-T3 antibody. Compound of biotinylated T3 and anti-T3 antibody is bridged to the microparticles by biotin. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction vessel. It is catalyzed by anti-T3 antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of free T3 present in the sample is inversely proportional to relative light units (RLUs) generated during the reaction. The Free T3 concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with streptavidin in HEPES buffer with preservative.
Rb	Monoclonal anti-T3 antibody (mouse)-alkaline phosphatase conjugate in MES buffer with preservative.
Rc	Biotinylated T3 in PBS buffer with preservative.
Rd	MES Buffer with preservative.

The position of each reagent component is shown in the figure below (front view on the left and top view on the right):

Storage and Stability

The unopened FT3 (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The FT3 (CLIA) reagent kit can be stored onboard and used for a maximum of 56 days after opening 2-8°C.

Reagent Preparation

Ra: Ready to use

Rb: Ready to use

Rc: Ready to use

Rd: Ready to use

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Cat. No. FT3211: Mindray Free T3 Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat. No. TFL321/TFL322/TFL323/TFL324: Mindray Thyroid Function Multi Control (L), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.

Cat. No. TFH321/TFH322/TFH323/TFH324: Mindray Thyroid Function Multi Control (H), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.

Cat.No.WB411: Mindray Wash Buffer, 1×10 L.

Cat. No. CS511/CS512: Mindray Substrate Solution, 4×115 mL/4×75 mL.

Mindray Reaction Vessels

Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Human serum is recommended for this assay.

Centrifuge the specimens after clot formation is complete. Transfer the supernatants into tubes for storage or test within two hours after centrifugation.

Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 72 hours, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles.

Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the FT3 (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to

use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 30 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

Calibration

CL-series FT3 (CLIA) has been standardized against a commercial FT3 test (CLIA).

The specific information of master calibration curve of FT3 (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any FT3 test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Mindray Thyroid Function Multi Control (L) and Thyroid Function Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of pg/mL.

Conversion factors:

pg/mL x 1.536 = pmol/L
pmol/L x 0.651 = pg/mL

Expected values

An extensive study on a cohort of 350 healthy individuals (175 males and 175 females) has determined the reference range of CL-series FT3 assay.

Category	N	Central 95% Range
Male	175	2.4-4.2 pg/mL
Female	175	2.4-3.9 pg/mL
Total	350	2.3-4.0 pg/mL

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

The upper limit of this assay is 30 pg/mL. A specimen with a free T3 concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >30 pg/mL.

The concentration of free T3 in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Performance Characteristics

Analytical Sensitivity/Limit of Detection

The FT3 (CLIA) reagent kit has an analytical sensitivity of ≤ 0.88 pg/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the FT3 concentration at two standard deviations below the mean RLU from 20 measurements of an analyte-free sample.

Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of FT3 (CLIA) reagent kit is 0.88-30 pg/mL.

Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 20 mg/dL and triglycerides up to 2000 mg/dL will not interfere with the CL-series FT3 assay. These substances show less than 10% interferences at indicated concentration.

In vitro tests were performed on 3 commonly used pharmaceuticals. These compounds

P/N: 046-003245-00 (7.0)

showed less than 10% interference in the CL-series FT3 assay at the levels indicated below.

Test Compound	Interference Substance Concentration
Methimazol	0.4 mg/dL
Phenylbutazone	1 mg/dL
Thiouracil	1 mg/dL

Mindray Free T3 Calibrator C0 were supplemented with other potential cross-reactants at specific levels indicated in the table below. No obvious cross reactivity was observed. The results are stated in the table below.

Substance	Cross-reactant Concentration	Cross Reactivity
Diiodo-L-tyrosine	76.0 pg/mL	0.00%
3,5 diiodothyronine	3.83 pg/mL	0.00%
Monoiodotyrosine	72.0 pg/mL	0.00%
Reverse T3	450.0 pg/mL	0.00%
Tetraiodothyroacetic acid	200,000 pg/mL	0.00%
L-thyroxine	50000 pg/mL	0.00%
3-3'5-Tri-iodothyroacetic Acid(TRIAC)	5000 pg/mL	0.29%

Accuracy

Two trueness controls with traceable and defined values were used to verify the accuracy of this assay. The results showed that the relative deviation was less than $\pm 10\%$. The results are listed in the following table.

Sample	Measured FT3 Value (pg/mL)	Defined FT3 Value (pg/mL)	Relative Deviation
Level 1	3.86	3.77	2.39%
Level 2	14.07	14.63	-3.83%

Precision

The CL-series FT3 assay is designed to have a precision of $\leq 10\%$ (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2. Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean FT3 (pg/mL)	Within-run CV	Between-run CV	Within-Device CV
1	4.62	2.65%	2.41%	3.31%
2	14.78	1.89%	2.16%	2.62%

Method Comparison


The Mindray CL-Series FT3 Assay was compared to a commercially available diagnostic kit in a correlation study with about 465 specimens. The statistical data obtained by Deming computing mode are shown in the table below.

Concentration Range (pg/mL)	Slope	Intercept	Correlation Coefficient
0.88-30	1.0124	0.4625	0.991

Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

Graphical Symbols

References

- B.S.Pan, Q.S.Fan, Q.Shen, J.P.Yu, M.Guan, Clinical application guide of protein detection in laboratories. Edition i. Shanghai Scientific & Technical Publishers, 2008.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.
- Ekins RP, editor. Methods for the Measurement of Free Thyroid Hormones. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
- Robbins J, Rall JE. The Iodine-Containing Hormones. In: Hormones in Blood (3rd Ed.). London: Academic Press, 1979:1:632-667.
- Gornall AG, Luxton AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318.
- White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
- Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE, Utiger RD, editors. The Thyroid (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
- Wahner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-26582888

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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- B.S.Pan, Q.S.Fan, Q.Shen, J.P.Yu, M.Guan, Clinical application guide of protein detection in laboratories. Edition i. Shanghai Scientific & Technical Publishers, 2008.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.

- Ekins RP, editor. Methods for the Measurement of Free Thyroid Hormones. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
- Robbins J, Rall JE. The Iodine-Containing Hormones. B: Hormones in Blood (3rd Ed.). London: Academic Press, 1979:1:632-667.
- Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318.
- White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
- Ladenson PW. Diagnosis of Thyrotoxicosis. B: Braverman LE, Utiger RD, editors. The Thyroid (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
- Wahner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943.

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Catalogue number

Referencias

- B.S.Pan, Q.S.Fan, Q.Shen, J.P.Yu, M.Guan, Clinical application guide of protein detection in laboratories. Edición i. Shanghai Scientific & Technical Publishers,2008.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.
- Ekins RP, editor. Methods for the Measurement of Free Thyroid Hormones. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
- Robbins J, Rall JE. The Iodine-Containing Hormones. En: Hormones in Blood (3ª ed.). London: Academic Press, 1979:1:632-667.
- Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318.
- White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously

- Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
- Ladenson PW. Diagnosis of Thyrotoxicosis. En: Braverman LE, Utiger RD, editors. The Thyroid (6th Ed.). Filadelfia: JB Lippincott Co., 1991:880-886.
 - Wahner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943.

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Manufacturer - Fabbricante

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- 1. B.S.Pan, Q.S.Fan, Q.Shen, J.P.Yu, M.Guan, Clinical application guide of protein detection in laboratories. Edition i. Shanghai Scientific & Technical Publishers, 2008.
- 2. Larsen PR. Triiodothyronine: Review of

- Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.
- 3. Ekins RP, editor. Methods for the Measurement of Free Thyroid Hormones. Amsterdam: Excerpta Medica Foundation. 1979;72-92.
- 4. Robbins J, Rall JE. The Iodine-Containing Hormones. In: Hormones in Blood (3rd Ed.). London: Academic Press, 1979:1:632-667.
- 5. Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318.
- 6. White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- 7. Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
- 8. Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE, Utiger RD, editors. The Thyroid (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
- 9. Wahner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943

© 2013 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd.
Tutti i diritti riservati
Produttore: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Indirizzo: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, Cina.

Indirizzo e-mail: service@mindray.com
Sito Web: www.mindray.com
Tel.: +86-755-26582888
Fax: +86-755-26582680
Rappresentante CE: Shanghai International Holding Corp. GmbH (Europe)
Indirizzo: Eiffestraße 80, Amburgo 20537, Germania
Tel.: 0049-40-2513175
Fax: 0049-40-255726

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
Grafiksel Semboller

Referanslar

- B.S.Pan, Q.S.Fan, Q.Shen, J.P.Yu, M.Guan, Clinical application guide of protein detection in laboratories. Edition i. Shanghai Scientific & Technical Publishers, 2008.
- Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.
- Ekins RP, editor. Methods for the Measurement of Free Thyroid

Hormones. Amsterdam: Excerpta Medica Foundation. 1979;72-92.

- Robbins J, Rall JE. The Iodine-Containing Hormones. In: Hormones in Blood (3rd Ed.). London: Academic Press, 1979;1:632-667.
- Gornall, AG, Luxton, AW, Bhavnani, BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986, 305-318.
- White, GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences, 1987, 24: 315-362.
- Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
- Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE, Utiger RD, editors. The Thyroid (6th Ed.). Philadelphia: JB Lippincott Co., 1991:880-886.
- Wahner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 Çin

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Symboles graphiques

In vitro diagnostic medical device

Batch code

European Conformity

Authorized representative in the European Community

Use by

Consult instructions for use

Caution

Temperature limit

Manufacturer

Catalogue number

Références

1. B.S.Pan. Q.S.Fan. Q.Shen. J.P.Yu. M.Guan. Clinical application guide of protein detection in laboratories. Edition i. Shanghai Scientific & Technical Publishers, 2008.
2. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. Metabolism 1972;21:1073-1092.
3. Ekins RP. editor. Methods for the Measurement of Free Thyroid Hormones.

Amsterdam: Excerpta Medica Foundation. 1979;72-92.

4. Robbins J. Rall JE. The Iodine-Containing Hormones. In: Hormones in Blood (3rd Ed.). London: Academic Press. 1979;1:632-667.
5. Gornall. AG. Luxton. AW. Bhavnani. BR. Endocrine disorders. In Applied Biochemistry of Clinical Disorders. 1986. 305-318.
6. White. GH. Recent advances in routine thyroid function testing. CRC - Critical Reviews in Clinical Laboratory Sciences. 1987. 24: 315-362.
7. Hamburger JL. Evolution of Toxicity in Solitary Nontoxic Autonomously Functioning Thyroid Nodules. J Clin Endocrinol Metab 1980;50:1089-1093.
8. Ladenson PW. Diagnosis of Thyrotoxicosis. In: Braverman LE. Utiger RD. editors. The Thyroid (6th Ed.). Philadelphia: JB Lippincott Co.. 1991:880-886.
9. Wahnner HW. T3 Hyperthyroidism. Mayo Clin Proc 1972;47:938-943.

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© 2013 Shenzhen Mindray Bio-Medical Electronics Co. Ltd.
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co. Ltd.
Représentant EC : Shanghai International Holding Corp. GmbH (Europe)