Total β HCG

Total β Human Chorionic Gonadotropin (CLIA)

Order Information

Catalog No.	Package Size
HCG111	2×50 tests
HCG112	2×100 tests

Intended Use

The CL-series Total β HCG assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of total β human chorionic gonadotropin (Total β HCG) in human serum.

Summary

Human chorionic Gonadotropin (HCG) is a glycoprotein hormone, produced by the placenta, with a molecular weight of approximately 38,000 daltons. It consists of 2 subunits (alpha-and beta-chains) that are associated to the intact hormone. The alpha subunit (MW 15,000-20,000 daltons) is common to FSH, TSH, and LH, but the beta subunit differ, and confer immunological and biological specificity. The β-chains (MW 25,000-30,000 daltons) have greatly differing structures and are responsible for the respective specific hormonal functions.^{1,2}

HCG is produced in the placenta during pregnancy. In non-pregnant women, it can also be produced by tumors of the trophoblast, germ cell tumors with trophoblastic components and some non-trophoblastic tumors.³ Human chorionic gonadotropin consists of a number of isohormones with differing molecular size. The biological action of HCG serves to maintain the corpus luteum during pregnancy. It also influences steroid production.

The serum of pregnant women contains mainly intact HCG. HCG molecule appears between 2 and 8 days after fertilization in the embryo's culture medium. HCG can be detected in maternal serum 7 days after fertilization and 8 days following ovulation.^{4,5,6} HCG increases 3-fold between the day of detection and the next day, decreasing there after and reaching 1.6-fold between Days 6 and 7. The HCG concentration increases exponentially in the first trimester, doubling every 31-48 h and reaches greater amounts between 11 and 13 weeks. In the second trimester, at 20 weeks gestation, it diminishes to 80% and remains in this concentration until the term. The HCG is more negatively charged in the early than in the latter stage of gestation and this change in charge occurs around Week 13, when the placental production of estradiol and progesterone increases. Unusually low or rapidly declining levels may indicate an abnormal condition such as an ectopic pregnancy or impending spontaneous abortion.

Assay Principle

The Total β HCG assay is a two-site sandwich assay to determine the level of Total β HCG.

In the first step, sample, paramagnetic microparticle coated with monoclonal anti-β HCG antibody (mouse) are added into a reaction vessel. After incubation, β HCG present in the sample binds

to anti-β HCG antibody coated microparticle. Microparticle is magnetically captured while other unbound substances are removed by washing.


In the second step, the monoclonal anti-β HCG antibody (mouse)-alkaline phosphatase conjugate is added into the reaction vessel. After incubation, the HCG binds to anti-β HCG antibody coated microparticle and anti-β HCG antibody alkaline phosphatase-labeled conjugate to form a sandwich complex. Microparticle is magnetically captured while other unbound substances are removed by washing.

In the third step, the substrate solution is added to the reaction vessel. It is catalyzed by anti- β HCG antibody (mouse)-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of β HCG present in the sample is proportional to the relative light units (RLUs) generated during the reaction. The β HCG concentration can be determined via a calibration curve.

Reagent Components

Ra	Paramagnetic microparticles coated with monoclonal anti- β HCG antibody (mouse) in TRIS buffer with preservative.
Rb	Monoclonal anti-βHCG antibody (mouse)-alkaline phosphatase conjugate in PBS buffer with preservative.
Rc	TRIS buffer with preservative.

The position of each reagent component is shown in the figure below (front view on the left and top view on the right):

Storage and Stability

The unopened Total β HCG (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C.

The Total β HCG (CLIA) reagent kit can be stored onboard and used for a maximum of 28 days after opening at 2-8°C.

Reagent Preparation

Ra: Ready to use

Rb: Ready to use

Rc: Ready to use

Materials Required but not Provided

Mindray CL-series Chemiluminescence Immunoassay Analyzer.

Cat.No.HCG211: Mindray Total β HCG Calibrators, 1×2.0 mL for each of calibrator C0, C1 and C2.

Cat.No.RHL311: Mindray Reproductive Multi Control (L), 6×5.0 mL.

Cat.No.RHH311: Mindray Reproductive Multi

Control (H), 6×5.0 mL.

Cat.No.RHL312: Mindray Reproductive Multi Control (L), 12×5.0 mL.

Cat.No.RHH312: Mindray Reproductive Multi Control (H), 12×5.0 mL.

Cat.No.RHL313: Mindray Reproductive Multi Control (L), 1×5.0 mL.

Cat.No.RHH313: Mindray Reproductive Multi Control (H), 1×5.0 mL.

Cat.No.RHL314: Mindray Reproductive Multi Control (L), 3×5.0 mL.

Cat.No.RHH314: Mindray Reproductive Multi Control (H), 3×5.0 mL.

Cat.No.WB411: Mindray Wash Buffer, 1×10 L.

Cat.No.CS511: Mindray Substrate Solution, 4×115 mL.

Cat.No.CS512: Mindray Substrate Solution, 4×75 mL.

Mindray Reaction Vessels

Applicable Instrument

Mindray CL-series Chemiluminescence Immunoassay Analyzer

Specimen Collection and Preparation

Human serum is recommended for this assay.

Centrifuge the specimens after clot formation is complete. Transfer the supernatants into tubes for storage or test within two hours after centrifugation.

Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 72 hours, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles.

Assay Procedure

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well.

Before loading the Total β HCG (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 10 μL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the

assay to determine the minimum sample volume.

Calibration

CL-series Total β HCG (CLIA) has been standardized against the WHO International Standard 5th WHO IS Chorionic Gonadotropin NIBSC code: 07/364.

The specific information of master calibration curve of Total β HCG (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any β HCG test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

Quality Control

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Mindray Reproductive Multi Control (L) and Reproductive Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

Calculation

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of mIU/mL.

Conversion factors: mIU/mL x 1 = IU/L

Dilution

Samples with Total β HCG concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:40 (either automatically by the analyzer or manually). The concentration of the diluted sample must be > 10 mIU/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiply the result by the dilution factor when calculating the sample concentration.

Expected values

An extensive study on a cohort of 451 healthy individuals (137 males and 314 non-pregnant females) has determined the reference range of CL-series β HCG assay.

Reference Group	N	Central 95% Range
Males	137	<0,5-2,2 mIU/mL

Non-pregnant females	314	<0,5-3,0 mIU/mL
----------------------	-----	-----------------

Representative β HCG ranges during normal pregnancy are summarized below.

Approximate Gestational Age (weeks)	N	Central 95% Range
1-10	129	50-221796 mIU/mL
11-15	121	15205-227034 mIU/mL
16-24	122	5002-75445 mIU/mL
25-40	121	1501-69208 mIU/mL

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

Limitation

This kit is NOT intended to be used for the risk evaluation of trisomy 21.

The upper limit of this assay is 5000mIU/mL.

A specimen with a β HCG concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >5000 mIU/mL or diluting the samples with Mindray Sample Diluent.

The concentration of β HCG in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies. However, no obvious interference of HAMA was observed in this assay.

Performance Characteristics

Analytical Sensitivity/Limit of Detection

The Total β HCG (CLIA) reagent kit has an analytical sensitivity of ≤ 0.5 mIU/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the β HCG concentration at two standard deviations above the mean RLU from 20 measurements of an analyte-free sample.

Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of Total β HCG (CLIA) reagent kit is 0.5-5000 mIU/mL (or the upper limit is up to 200000 mIU/mL for 40-fold diluted samples).

Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 10 mg/dL, triglycerides up to 1400 mg/dL, and total protein up to 10.0 g/dL will not interfere with the CL-series β HCG assay. These substances show less than 10% interferences at indicated concentration.

No obvious interference was observed from rheumatoid factor up to 800 IU/mL or antinuclear

antibody up to 2000 U/L.

Mindray Total HCG Calibrator C0 was supplemented with follicle stimulating hormone (FSH), thyroid stimulating hormone (TSH), or luteinizing hormone (LH), at specific levels indicated in the table below. No obvious cross reactivity was observed as all the results were ≤ 3.0 mIU/mL. The results are stated in the table below.

Sub-stance	Cross-reactant Concentration	Reported β HCG (mIU/mL)	Acceptance Criteria
FSH	200 mIU/mL	0,14	Reported β HCG $\leq 3,0$ mIU/mL
TSH	100 μ IU/mL	0,25	
LH	200 mIU/mL	0,22	

Accuracy

A sample of WHO International Standard 5th WHO IS Chorionic Gonadotropin (NIBSC code: 07/364) with traceable and defined value was used to verify the accuracy of this assay. The results showed that the relative deviation was less than $\pm 10\%$. The results are listed in the following table.

Sample	Measured β HCG Value (mIU/mL)	Defined β HCG Value (mIU/mL)	Relative Deviation
WHO HCG	24,07	25,80	-6,66%

Precision

The CL-series Total β HCG assay is designed to have a precision of $\leq 10\%$ (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2. Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean β HCG (mIU/mL)	Within-run CV	Between-run CV	Within-Device CV
1	6,72	2,50%	1,38 %	2,88 %
2	191,01	1,14%	0,73 %	1,26 %

Linearity

A high concentration β HCG sample (approximately 5000 mIU/mL) was mixed with a low concentration sample (< 5 mIU/mL) at different ratios, generating a series of dilutions. The β HCG of each dilution was determined using the Mindray CL-Series β HCG Assay. Linearity was demonstrated in the range of 5.0 mIU/mL to 5000 mIU/mL. The correlation coefficient r is ≥ 0.9900 . The linearity data are summarized in the table below.

Concentration (mIU/mL)	1	2	3	4	5	6
Expected β HCG	0,4	1008,40	2016,76	3025,12	4033,47	5041,83
Measured β HCG	0,4	1181,67	2255,73	3485,17	4047,27	5041,83

Method Comparison

The Mindray CL-Series Total β HCG Assay


was compared to a commercially available diagnostic kit in a correlation study with about 447 specimens. The statistical data obtained by Deming computing mode are shown in the table below.

Concentration Range (mIU/mL)	Slope	Intercept	Correlation Coefficient
0,5-1000	1,04	3,57	0,9903

Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 28 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

Graphical Symbols

References

- Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Issacs N, Jackson A, Klonisch T, Laphorn A, Lund T et al. Immunochemical mapping of gonadotropins. Mol Cell Endocrinol 1996;125:33-43.
- Bidart JM, Birken S, Borger P, Krichevsky A. Immunochemical mapping of hCG and hCG-related molecules. Scand J Clin Lab Invest 1993;53:118-136.
- Borrelli PTA, Butler SA, Docherty SM, Staite EM, Borrelli AL, Iles RK. Human chorionic gonadotropin isoforms in the diagnosis of ectopic pregnancy. Clin Chem 2003;12:2045-2049.

- Butler SA, Iles RK. The free monomeric beta subunit of human chorionic gonadotropin (hCG beta) and the recently identified homodimeric beta-beta subunit both have autocrine growth effects. Tumour Biol 2004;25:18-23.
- Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997;43:2233-2243.
- Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. Journal of clinical Immunoassay, Fall, 1991; 14, No. 3: 196-199.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved.

Manufacturer: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Address: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

E-mail Address: service@mindray.com

Website: www.mindray.com

Tel: +86-755-26582888

Fax: +86-755-26582680

EC-Representative: Shanghai International Holding Corp. GmbH (Europe)

Address: Eiffestraße 80, Hamburg 20537, Germany

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Símbolos gráficos

In vitro diagnostic medical device

Batch code

European Conformity

Authorized representative in the European Community

Use by

Consult instructions for use

Caution

Temperature limit

Manufacturer

Catalogue number

Referências

1. Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Issacs N, Jackson A, Klonisch T, Lapthorn A, Lund T et al. Immunochemical mapping of gonadotropins. Mol Cell Endocrinol 1996;125:33-43.
2. Bidart JM, Birken S, Borger P, Krichevsky A. Immunochemical mapping of hCG and hCG-related molecules. Scand J Clin Lab Invest 1993;53:118-136.
3. Borrelli PTA, Butler SA, Docherty SM, Staite EM, Borrelli AL, Iles RK. Human chorionic gonadotropin isoforms in the diagnosis of ectopic pregnancy. Clin Chem 2003;12:2045-2049.
4. Butler SA, Iles RK. The free monomeric beta subunit of human chorionic gonadotropin (hCG beta) and the recently identified homodimeric beta-beta subunit both have autocrine growth effects. Tumour Biol 2004; 25:18-23.
5. Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997;43:2233-2243.
6. Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. Journal of clinical Immunoassay, Fall, 1991; 14, No. 3: 196-199.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Fabricante: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Endereço: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, P.R. China.

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Catalogue number

Referencias

- Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Issacs N, Jackson A, Klonisch T, Laphorn A, Lund T et al. Immunochemical mapping of gonadotropins. Mol Cell Endocrinol 1996;125:33-43.
- Bidart JM, Birken S, Borger P, Krichevsky A. Immunochemical mapping of hCG and hCG-related molecules. Scand J Clin Lab Invest 1993;53:118-136.
- Borrelli PTA, Butler SA, Docherty SM, Staite EM, Borrelli AL, Iles RK. Human chorionic gonadotropin isoforms in the diagnosis of ectopic pregnancy. Clin Chem 2003;12:2045-2049.
- Butler SA, Iles RK. The free monomeric beta subunit of human chorionic gonadotropin (hCG beta) and the recently identified homodimeric beta-beta subunit both have autocrine growth effects. Tumour Biol 2004;25:18-23.
- Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997;43:2233-2243.
- Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. Journal of clinical Immunoassay, Fall, 1991; 14, No. 3: 196-199.

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Manufacturer - Fabbricante

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- Berger P, Bidart JM, Delves PS, Dirnhofer S, Hoermann R, Issacs N, Jackson A, Klonisch T, Laphorn A, Lund T et al. Immunochemical mapping of gonadotropins. Mol Cell Endocrinol 1996;125:33-43.
- Bidart JM, Birken S, Borger P, Krichevsky A. Immunochemical mapping of hCG and hCG-related molecules. Scand J Clin Lab Invest 1993;53:118-136.
- Borrelli PTA, Butler SA, Docherty SM, Staite EM, Borrelli AL, Iles RK. Human chorionic gonadotropin isoforms in the diagnosis of ectopic pregnancy. Clin Chem 2003;12:2045-2049.
- Butler SA, Iles RK. The free monomeric beta subunit of human chorionic gonadotropin (hCG beta) and the recently identified homodimeric beta-beta subunit both have autocrine growth effects. Tumour Biol 2004;25:18-23.
- Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997;43:2233-2243.
- Sokolove PJ, Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. Journal of Clinical Pathology, 1991, 14: 3: 196-199.

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Indirizzo: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, Cina.

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Catalogue number

Références

1. Berger P. Bidart JM. Delves PS. Dirnhofer S. Hoermann R. Issacs N.Jackson A. Klonisch T. Laphorn A. Lund T et al. Immunochemical mapping of gonadotropins.Mol Cell Endocrinol 1996;125:33-43.
2. Bidart JM. Birken S. Borger P. Krichevsky A. Immunochemical mapping of hCG and hCG-related molecules. Scand J Clin Lab Invest 1993;53:118-136.
3. Borrelli PTA. Butler SA. Docherty SM. Staite EM. Borrelli AL. Iles RK.Human chorionic gonadotropin isoforms in the diagnosis of ectopic pregnancy. Clin Chem 2003;12:2045-2049.
4. Butler SA. Iles RK. The free monomeric beta subunit of human chorionic gonadotropin (hCG beta) and the recently identified homodimeric beta-beta subunit both have autocrine growth effects. Tumour Biol 2004;25:18-23.
5. Cole LA. Immunoassay of human chorionic gonadotropin, its free subunits, and metabolites. Clin Chem 1997;43:2233-2243.
6. Sokolove PJ. Faix JD. Agreement of intact and beta chain-specific HCG assays in abnormal pregnancy. Journal of clinical Immunoassay. Fall. 1991; 14. No. 3: 196-199.

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© 2013 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd.
Fabricant : Shenzhen Mindray Bio-Medical Electronics Co., Ltd.
Représentant EC : Shanghai International Holding Corp. GmbH (Europe)