**AFP****Alpha-fetoprotein (CLIA)****Order information**

Catalog No.	Package Size
AFP111	2×50 tests
AFP112	2×100 tests

**Intended Use**

The CL-series AFP assay is a Chemiluminescent Immunoassay (CLIA) for the quantitative determination of alpha-fetoprotein (AFP) in human serum or plasma.

**Summary**

Alpha-fetoprotein (AFP) is an embryo-specific alpha-globulin and a single polypeptide chain glycoprotein with a molecular weight of approximately 70,000 Daltons. It is the main protein in the fetal circulation, as well as albumin. Synthesis of AFP occurs primarily in the yolk sac, liver and gastro-intestinal tract of the fetus.<sup>1</sup> It is secreted into fetal serum and can be detected in fetal circulation as early as 30 days after gestation. It reaches a peak at about 13 weeks gestation and gradually declining thereafter.<sup>2, 3</sup> Elevated serum AFP levels subsequently reappear during pregnancy and in conjunction with several malignant diseases.<sup>2, 4</sup>

Elevated serum AFP levels can be detected in about 80% patients of hepatic carcinoma and 50% patients of germ cell tumors. It is also can be detected in other patients such as patients of gastro-intestinal tumor. Synthesis of AFP is obviously elevated in hepatic carcinoma even though the corpus carcinosus is very small. So, hepatic carcinoma can be detected as early as possible by AFP detection, in conjunction with some other tests.<sup>5, 6</sup> Serum AFP has also been found to be elevated during pregnancy, in diseases such as acute viral hepatitis, chronic active hepatitis and cirrhosis, but not as high as hepatic carcinoma. The level of serum AFP in patients of cirrhosis is generally between 25 and 200 ng/mL, and gradually declines with the growing well of patients' condition with 2 months. The elevated level of serum AFP generally accompanies an elevate level of transaminase. If the level of serum AFP is elevated while the level of transaminase is declined, or the level of serum AFP is more than 500 ng/mL, it is highly possible to be hepatic carcinoma.<sup>6,7,8</sup> AFP testing is not recommended as a screening procedure to detect cancer in the general population.

**Assay Principle**

The CL-series AFP assay is a two-site sandwich assay to determine the level of AFP.

In the first step, sample, paramagnetic microparticles coated with mouse monoclonal anti-AFP antibody and mouse monoclonal anti-AFP antibody alkaline phosphatase labeled conjugate are added into a reaction vessel. After incubation, AFP present in the sample binds to the mouse monoclonal anti-AFP antibody coated microparticles. Meanwhile, the mouse monoclonal anti-AFP


antibody alkaline phosphatase labeled conjugate binds to another reaction site of the AFP to form a sandwich complex. Microparticles are magnetically captured while other unbound substances are removed by washing.

In the second step, the substrate solution is added to the reaction vessel. It is catalyzed by mouse monoclonal anti-AFP antibody-alkaline phosphatase conjugate in the immunocomplex retained on the microparticle. The resulting chemiluminescent reaction is measured as relative light units (RLUs) by a photomultiplier built into the system. The amount of AFP present in the sample is proportional to relative light units (RLUs) generated during the reaction. The AFP concentration can be determined via a calibration curve.

**Reagent Components**

Ra	Paramagnetic microparticles coated with monoclonal anti-AFP antibody (mouse) in TRIS buffer with preservative.
Rb	Monoclonal anti-AFP antibody (mouse)-alkaline phosphatase conjugate in PBS buffer with preservative.
Rc	TRIS buffer with preservative.

The position of each reagent component is shown in the figure below (front view on the left and top view on the right):

**Storage and Stability**

The unopened AFP (CLIA) reagent kit is stable up to the stated expiration date when stored at 2-8°C. The AFP (CLIA) reagent kit can be stored onboard and used for a maximum of 56 days after opening at 2-8°C.

**Reagent Preparation**

Ra: Ready to use  
Rb: Ready to use  
Rc: Ready to use

**Materials Required But Not Provided**

Mindray CL-series Chemiluminescence Immunoassay Analyzer  
Catalog No. AFP211: Mindray AFP Calibrators. 1×2.0 mL for each of calibrator C0, C1 and C2.  
Cat.No.TML311/TML312/TML313/TML314: Mindray Tumor Marker Multi Control (L), 6×5.0 mL/12×5.0 mL/1×5.0 mL/3×5.0 mL.  
Cat.No.TMH311/TMH312/TMH313/TMH314: Mindray Tumor Marker Multi Control (H), 6×5.0 mL/12×5.0 mL/1×5.0 mL/3×5.0 mL.  
Cat.No.TML321/TML322/TML323/TML324: Mindray Tumor Marker Multi Control (L), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.  
Cat.No.TMH321/TMH322/TMH323/TMH324: Mindray Tumor Marker Multi Control (H), 1×5.0 mL/3×5.0 mL/6×5.0 mL/12×5.0 mL.  
Cat.No.WB411: Mindray Wash Buffer, 1×10 L.

Cat.No.CS511/CS512: Mindray Substrate Solution, 4×115 mL/4×75 mL.

Mindray Reaction Vessels

**Applicable Instrument**

Mindray CL-series Chemiluminescence Immunoassay Analyzer

**Specimen Collection and Preparation**

Human serum or plasma collected in EDTA, sodium heparin and lithium heparin are recommended for this assay.

Centrifuge the specimens after clot formation is complete. Transfer the supernatants into tubes for storage or test within two hours after centrifugation. Specimens should be tested as soon as possible after sample collection. If testing is not completed within 8 hours, specimens should be tightly capped and refrigerated at 2-8°C. If testing will be delayed for more than 72 hours, specimens should be frozen at -20°C or below.

Avoid repeated freeze and thaw cycles.

**Assay Procedure**

For optimal performance of this assay, operators should read the related system operation manual carefully, to get sufficient information such as operation instructions, sample preservation and management, safety precaution, and maintenance. Prepare all required materials for the assay as well. Before loading the AFP (CLIA) reagent kit on the machine for the first time, unopened reagent bottle should be inverted gently for at least 30 times to resuspend the microparticles that have settled during shipment or storage. Visually inspect the bottle to ensure the microparticles have been resuspended. If the microparticles remain adhered to the bottle, continue inverting until the microparticles have been completely resuspended. If the microparticles cannot be resuspended, it is recommended not to use this bottle of reagent. Contact Mindray Customer Service for help. Do not invert opened reagent bottle.

This assay requires 15 µL of sample for a single test. This volume does not include the dead volume of the sample container. Additional volume is required when performing additional tests from the same sample. Operators should refer to the system operation manual and specific requirement of the assay to determine the minimum sample volume.

**Calibration**

CL-series AFP (CLIA) has been standardized against the WHO International Standard Alphafoetoprotein (NIBSC code: AFP).

The specific information of master calibration curve of AFP (CLIA) reagent kit is stored in the two-dimensional barcode attached in the reagent pack. It's used together with calibrators for the calibration of the specific reagent lot. When performing the calibration, scan the information of master calibration curve from the barcode into the system first, and then use the calibrators at three levels. Valid calibration curve is required before any AFP test. Recalibration is recommended every 4 weeks, or when a new reagent lot is used, or the

quality controls are out of specified range. For detailed instruction of calibration, refer to the system operation manual.

**Quality control**

It is recommended that quality controls should be run once every 24 hours if the tests are in use, or after every calibration. The quality control frequency should be adapted to each laboratory's individual requirements. The recommended two levels of quality controls for this assay are Mindray Tumor Marker Multi Control (L) and Tumor Marker Multi Control (H).

Quality control results should be within the acceptable ranges. If a control is out of its specified range, the associated test results are invalid and the samples must be retested. Recalibration may be required. Examine the assay system referring to the system operation manual. If the quality control results are still out of the specified range, please contact Mindray Customer Service for help.

**Calculation**

The analyzer automatically calculates the analyte concentration of each sample on the master calibration curve read from the barcode, and a 4-Parameter Logistic Curve Fitting (4PLC) with the relative light units (RLUs) generated from three level calibrators of defined concentration values. The results are shown in the unit of ng/mL.  
Conversion factors: ng/mL × 0.83 = IU/mL  
IU/mL × 1.21 = ng/mL

**Dilution**

Samples with AFP concentrations above the upper limit can be diluted with Mindray Sample Diluent. The recommended dilution is 1:40 (either automatically by the analyzer or manually). The concentration of the diluted sample must be > 10 ng/mL. After manual dilution, multiply the result by the dilution factor. After automated dilution by the analyzers, the system automatically multiply the result by the dilution factor when calculating the sample concentration.

**Expected values**

An extended study on a cohort of 736 healthy individuals (345 males and 391 females) has determined the reference range of CL-series AFP assay.

Category	N	97.5% Percentile
Male	345	7.2 ng/mL
Female	391	7.4 ng/mL
Total	736	7.0 ng/mL

Due to the variation in geography, race, sex, and age, it is highly recommended that each laboratory should establish its own reference range.

**Limitation**

This kit is NOT intended to be used for the risk evaluation of trisomy 21.

The upper limit of this assay is 1210 ng/mL. A specimen with a AFP concentration lower than the upper limit can be quantitatively determined, while specimen with a concentration higher than the upper limit will be reported as >1210 ng/mL or diluting the samples with Mindray Sample Diluent.

The concentration of AFP in a given specimen, determined with assays from different manufacturers, can vary due to differences in assay methods, calibration, and reagent specificity. The assay results should be used in conjunction with other data, such as symptoms, results of other tests, clinical history, etc.

Specimen from individuals who have been exposed to mouse monoclonal antibodies may contain human anti-mouse antibodies (HAMA). Such specimens may show either falsely elevated or depressed values with assay kits employing mouse monoclonal antibodies. However, no obvious interference of HAMA was observed in this assay.

### Performance Characteristics

#### Analytical Sensitivity/Limit of Detection

The AFP (CLIA) reagent kit has an analytical sensitivity of  $\leq 0.5$  ng/mL. Analytical sensitivity is defined as the lowest concentration of analyte that can be differentiated from a sample that contains no analyte. It is defined as the AFP concentration at two standard deviations above the mean RLU from 20 measurements of an analyte-free sample.

#### Reportable Range

Reportable range is defined by the analytical sensitivity and the upper limit of the master calibration curve. The reportable range of AFP (CLIA) reagent kit is 0.5-1210 ng/mL (or the upper limit is up to 48400 ng/mL for 40-fold diluted samples).

#### Specificity

Hemoglobin up to 500 mg/dL, bilirubin up to 20 mg/dL, triglycerides up to 520 mg/dL, and total protein up to 10.0 g/dL will not interfere with the CL-series AFP assay. These substances show less than 10% interferences at indicated concentration. No obvious interference was observed from rheumatoid factor up to 1500 IU/mL or antinuclear antibody up to 4000 U/L.

No significant cross-reactivity was observed when human chorionic gonadotropin (HCG),  $\alpha$ -1 acid glycoprotein, transferrin,  $\alpha$ -1 antitrypsin and prolactin (PRL) were added to the AFP Calibrator C0, at specific levels indicated in the table below. The results showed no more than 2.0 ng/mL in the CL-series AFP assay at the levels indicated below.

Substrates	Cross-reactant Concentration
HCG	1000 mIU/mL
$\alpha$ -1 acid glycoprotein	2 mg/mL
Transferrin	400 mg/ dL
$\alpha$ -1 antitrypsin	500 mg/dL
PRL	50,000 mIU/L

Interference substance tests were performed on 6 commonly used pharmaceuticals. These compounds showed less than 10% interference in the CL-series AFP assay at the levels indicated below.

Test Compound	Interference Substance Concentration
Acetaminophen	200 $\mu$ g/mL

Acetylsalicylic acid	0.50 mg/mL
Cisplatin	100 $\mu$ g/mL
Ascorbic acid	1000 $\mu$ g/mL
Cyclophosphamide	25 mg/ dL
Vincristine	70 mg/ dL

Mindray AFP Calibrator C0 was supplemented with other tumor markers, such as Prostate-Specific Antigen (PSA), cancer antigen 125 (CA125), cancer antigen 15-3 (CA15-3), carbohydrate antigen 19-9 (CA19-9), carcinoembryonic antigen (CEA) and ferritin (FERR) at specific levels indicated in the table below. No obvious cross reactivity was observed as all the results were  $\leq 5.0$  ng/mL. The results are stated in the table below.

Tumor Marker	Cross-reactant Concentration	Reported AFP (ng/mL)	Acceptance Criteria
PSA	100 ng/mL	0.00	Reported AFP $\leq 5.0$ ng/mL
CA125	1000 U/mL	0.01	
CA15-3	100 U/mL	0.01	
CA19-9	1000 U/mL	0.00	
CEA	1000 ng/mL	0.00	
FERR	1000 ng/mL	4.61	

#### High Dose Hook

For the CL-series AFP assay, no high dose hook effect was observed when samples containing up to approximately 120,000 ng/mL of AFP were assayed.

#### Accuracy

A sample of WHO International Standard Alphafoetoprotein (NIBSC code: AFP) with traceable and defined value was used to verify the accuracy of this assay. The results showed that the relative deviation was less than  $\pm 10\%$ . The results are listed in the following table.

Sample	Measured AFP Value (ng/mL)	Defined AFP Value (ng/mL)	Relative Deviation
WHO AFP	112.64	116.83	-3.59%

#### Precision

The CL-series AFP assay is designed to have a precision of  $\leq 10\%$  (within-device CV). Precision was determined by following National Committee for Clinical Laboratory Standards (NCCLS) Protocol EP5-A2.9 Two levels of quality controls were tested in duplicate in two separate runs per day, for a total of 20 days, using a single lot of reagents and a single calibration curve. The precision data are summarized in the table below.

Sample	Mean AFP (ng/mL)	Within-run CV	Between-run CV	Within-Device CV
1	10.32	2.92%	3.17 %	4.57 %
2	104.22	2.99%	2.52 %	3.98 %

#### Linearity

A high concentration AFP sample (approximately 400 ng/mL) was mixed with a low concentration

sample ( $< 1.0$  ng/mL) at different ratios, generating a series of dilutions. The AFP of each dilution was determined using the Mindray CL-Series AFP Assay. Linearity was demonstrated in the range of 1.0 ng/mL to 400 ng/mL. The correlation coefficient  $r$  is  $\geq 0.9900$ . The linearity data are summarized in the table below.

Concentration (ng/mL)	1	2	3	4	5	6
Expected AFP	0.14	104.75	209.35	313.96	418.56	523.16
Measured AFP	0.14	103.14	208.47	313.14	415.70	523.16

#### Method Comparison


The Mindray CL-Series AFP Assay was compared to a commercially available diagnostic kit in a correlation study with about 1287 specimens. The statistical data obtained by Deming computing mode are shown in the table below.

Concentration Range (ng/mL)	Slope	Intercept	Correlation Coefficient
0.13-1476.15	0.949	0.115	0.996

#### Warnings and Precautions

- For in vitro diagnostic use only.
- Follow all the rules in handling laboratory reagents and take necessary safety precautions.
- Due to the differences in methodology and antibody specificity, test results of the same sample may be different when using reagent kits from different manufacturers on Mindray system, or using Mindray reagent kits on other systems.
- Do not use reagent kits beyond the expiration date.
- Do not use reagents mixed from different reagent lots.
- Always keep the reagent pack in the upright position to ensure no microparticle has been lost prior to use.
- Reagent pack opened for more than 56 days is not recommended for use.
- Reliability of assay results cannot be guaranteed if the instructions in this package insert are not followed.
- All the specimen and reaction wastes should be considered potentially biohazard. The handling of specimens and reaction wastes should be in accordance with the local regulations and guidelines.
- The Material Safety Data Sheet (MSDS) is available upon request.

#### Graphical symbols

#### References

- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of

human  $\alpha$ -fetoprotein. Int. J. Cancer 1971; 7: 218.

- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastroschisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
- Gitlin D, Perricelli A, Gitlin GM. Synthesis of  $\alpha$ -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. In: Sull SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
- Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
- Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
- Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of  $\alpha$ -fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

All rights Reserved

**Manufacturer:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Address:** Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 P.R.China

**E-mail Address:** service@mindray.com

**Website:** www.mindray.com

**Tel:** +86-755-26582888

**Fax:** +86-755-26582680

**EC-Representative:** Shanghai International Holding Corp. GmbH (Europe)

**Address:** Eiffestraße 80, Hamburg 20537, Germany

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- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human a-fetoprotein. Int. J. Cancer 1971; 7: 218.
- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastroschisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
- Gitlin D, Perricelli A, Gitlin GM. Synthesis of a-fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1977; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. In: Sell SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
- 6 Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
- Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
- Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of a-fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

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
#### Símbolos gráficos

#### Referências

1. Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human  $\alpha$ -fetoprotein. Int. J. Cancer 1971; 7: 218.
2. Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastroschisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
3. Gitlin D, Perricelli A, Gitlin GM. Synthesis of  $\alpha$ -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
4. Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
5. Taketa K. Alpha-Fetoprotein in the 1990s. Em: Sell SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
6. Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
7. Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
8. Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of  $\alpha$ -fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
9. NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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**Fabricante:** Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

**Endereço:** Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, P.R. China.

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
##### Símbolos gráficos

##### Referencias

- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human a-fetoprotein. Int. J. Cancer 1971; 7: 218.
- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastroschisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
- Gitlin D, Perricelli A, Gitlin GM. Synthesis of a-fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. En: Sell SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
- Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
- Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
- Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of a-fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

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**Materiali necessari ma non forniti**

Analizzatore per immunodosaggio chemiluminescente Mindray serie CL  
Mindray AFP Calibrators, N. cat. AFP211: 1x2,0 ml per ciascun calibratore C0, C1 e C2.  
Mindray Tumor Marker Multi Control (L) N. cat. TML311/TML312/TML313/TML314: 6x5,0 ml/12x5,0 ml/1x5,0 ml/3x5,0 ml.  
Mindray Tumor Marker Multi Control (H) N. cat. TMH311/TMH312/TMH313/TMH314: 6x5,0 ml/12x5,0 ml/1x5,0 ml/3x5,0 ml.  
Mindray Tumor Marker Multi Control (L) N. cat. TML321/TML322/TML323/TML324: 1x5,0 ml/3x5,0 ml/6x5,0 ml/12x5,0 ml.  
Mindray Tumor Marker Multi Control (H) N. cat. TMH321/TMH322/TMH323/TMH324: 1x5,0 ml/3x5,0 ml/6x5,0 ml/12x5,0 ml.  
Tampone di lavaggio Mindray, 1x10 L, N. cat. WB411  
Soluzione substrato Mindray, 4x115 mL N. cat. CS511.  
Soluzione substrato Mindray, 4x75 mL N. cat. CS512.

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**Use By** – Scadenza

**Consult Instructions for Use** – Leggere

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**Temperature Limit** – Temperatura di

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**Manufacturer** – Fabbricante

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- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human  $\alpha$ -fetoprotein. Int. J. Cancer 1971; 7: 218.
- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastroschisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
- Gitlin D, Perricelli A, Gitlin GM. Synthesis of  $\alpha$ -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. In: Sell SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
- 6 Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
- Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
- Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of  $\alpha$ -fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

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**Indirizzo:** Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, Cina.

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
#### Grafiksel Semboller

#### Referanslar

- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human  $\alpha$ -fetoprotein. Int. J. Cancer 1971; 7: 218.
- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastrochisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.

- Gitlin D, Perricelli A, Gitlin GM. Synthesis of  $\alpha$ -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. In: Sell SS. Serological cancer markers. Humana Press 1992: 31-46, ISBN: 0-89603-209-4.
- Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human alpha-fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
- Perlin E, Engeler JE, Edson M, Karp D, McIntire KR, Waldmann TA. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. Cancer 1976; 37: 215-9.
- Endo Y, Kanai K, Oda T, Mitamura K, Iino S, Suzuki H. Clinical significance of  $\alpha$ -fetoprotein in hepatitis and liver cirrhosis. Ann NY Acad Sci 1975; 259: 234.
- NCCLS. EP5-A2 Evaluation of precision performance of quantitative measurement method; approved guideline-second edition. 2004.

© 2013 Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

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Üretici: Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Adres: Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, Shenzhen, 518057 Çin

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
##### Symboles graphiques

##### Références

- Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human  $\alpha$ -fetoprotein. Int. J. Cancer 1971; 7: 218.
- Palomaki GE, Hill LE, Knight GJ, et al. Second Trimester Maternal Serum Alpha-Fetoprotein Levels in Pregnancies Associated with Gastrochisis and Omphalocele. Obstet Gynecol 1988; 71:906-9.
- Gitlin D, Perricelli A, Gitlin GM. Synthesis of  $\alpha$ -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. Cancer Res 1972; 32: 979.
- Nelson LH, Bensen J, Burton BK. Outcomes in Patients with Unusually High Maternal Serum Alpha-fetoprotein Levels. Am J Obstet Gynecol 1987; 157(3):572-5.
- Taketa K. Alpha-Fetoprotein in the 1990s. In: Sell SS. Serological cancer markers. Humana Press 1992: 31-46. ISBN: 0-89603-209-4.
- Silver HKB, Gold P, Feder S, Freedman SO, Shuster J. Radioimmunoassay for human

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$\alpha$ -fetoprotein. Proc Natl Acad Sci U.S.A. 1973; 70: 526-30.
© 2013 Shenzhen Mindray Bio-Medical Electronics Co. Ltd.
**Fabricant :** Shenzhen Mindray Bio-Medical Electronics Co. Ltd.
**Représentant EC :** Shanghai International Holding Corp. GmbH (Europe)